

MDA5

BLUF: lineblot is highly prone to FP and FN results. IP is best overall, but ELISA is a faster, acceptable non-IP alternative.

Modality ranked by highest to lower *specificity*:

- IP > ELISA » multiplex bead assays » lineblot

Modality ranked by highest to lower *sensitivity*:

- ELISA >= IP » multiplex bead assays » lineblot

From:

<https://ewrobbins.com/> - ewrobbins.com

Permanent link:

https://ewrobbins.com/doku.php?id=ctd_testing&rev=1790178838

Last update: **2026/09/23 15:53**

