

SARD-ILD

ELISA

JO1, TIF1G, MDA5, NXP2, RO52, U1RNP

- item :: LabCorp, test number 520085. MyoMarker® 3 Plus Profile (RDL)
- instruction :: collect \geq 5 mL of blood in Red-top tube or gel-barrier tube

From:

<https://ewrobbins.com/> - **ewrobbins.com**

Permanent link:

https://ewrobbins.com/doku.php?id=lab_testing&rev=1778680120

Last update: **2026/05/13 13:48**

